

The low back W1i Tillett seat is supplied either as a single skin, or with a self supporting floor mountable back frame. If bought as a single skin shell it must have brackets made in order to fit it to various vehicles. These seats can also be ordered with a carbon or GRP front surface and with or without covering. A Dinamica suede panel set is also available. Three seat belt harness holes are included as standard. When ordered with the back frame you can select a "side mount" or "bottom mount" version. When side mounts are selected the TBFIA, STBi, STBO, VTBi, VTBO can be used. **Warning** Always use a load spreading spacer of at least 40 mm diameter when mounting. With bottom mounts do not test by sitting the seats on long bolts as you will potentially force the internal plate away from the laminate.

W1i-44 Carbon GRP with back frame



W1i-44 GRP shell only



W1i-40 GRP with back frame



Seat above and below is shown with the optional back frame.



The gap between the lowest point of the dome and the bolting surface is 38mm. With a bottom mount seat this gap must be filled with a spacer. Either the STRi / PTRi runners, and a TK5 spacer kit, or it can be set on aluminium or steel box section.

This model is not FIA approved.